

## GUARNIFLON S.p.A. - PATI Division

Product code : 0TQHST9900 - 100 mic

Direction : TD

Preload : 0,05 MPa


Rate of grip separation : 500 mm/min

Initial grip separation : 50,00 mm

### Risultati del test:

Sample no.	Thickness $\mu\text{m}$	$S_0$ $\text{mm}^2$	$E_t$ 0.1% MPa	Stress 10% MPa	Tensile strength at break MPa	Tensile strain at break %
1	101	1,01	767	18,1	48,4	340
2	101	1,01	757	18,1	47,8	340
3	102	1,02	752	17,8	47,7	340
4	105	1,05	819	19,4	52,3	350
5	101	1,01	775	18,3	50,3	350

### Grafico della serie:



### Statistiche:

Serie	Thickness $\mu\text{m}$	$S_0$ $\text{mm}^2$	$E_t$ 0.1% MPa	Stress 10% MPa	Tensile strength at break MPa	Tensile strain at break %
n = 5						
$\bar{x}$	102	1,02	774	18,3	49,3	340