

## GUARNIFLON S.p.A. - PATI Division

Product code : OTQHST9900 - 200 mic

Direction : MD

Preload : 0,05 MPa


Rate of grip separation : 500 mm/min

Initial grip separation : 50,00 mm

### Risultati del test:

Sample no.	Thickness $\mu\text{m}$	$S_0$ $\text{mm}^2$	$E_t$ 0.1% MPa	Stress 10% MPa	Tensile strength at break MPa	Tensile strain at break %
1	199	1,99	898	21,5	58,9	380
2	195	1,95	923	22,4	50,6	330
3	199	1,99	957	24,0	54,3	340
4	192	1,92	944	23,0	53,2	340

### Grafico della serie:



### Statistiche:

Serie	Thickness $\mu\text{m}$	$S_0$ $\text{mm}^2$	$E_t$ 0.1% MPa	Stress 10% MPa	Tensile strength at break MPa	Tensile strain at break %
n = 4						
$\bar{x}$	196,3	1,96	930	22,7	54,2	350