

GUARNIFLON S.p.A. - PATI Division

Product code : OTQHST9900 - 250 mic

Direction : MD

Preload : 0,05 MPa


Rate of grip separation : 500 mm/min

Initial grip separation : 50,00 mm

Risultati del test:

Sample no.	Thickness μm	S_0 mm^2	E_t 0.1% MPa	Stress 10% MPa	Tensile strength at break MPa	Tensile strain at break %
1	252	2,52	919	24,1	51,0	370
2	258	2,58	886	23,7	50,5	370
3	251	2,51	875	23,6	50,9	370
4	255	2,55	870	23,2	51,4	370
5	260	2,60	870	23,9	50,6	380

Grafico della serie:



Statistiche:

Serie	Thickness μm	S_0 mm^2	E_t 0.1% MPa	Stress 10% MPa	Tensile strength at break MPa	Tensile strain at break %
n = 5						
\bar{x}	255,2	2,55	884	23,7	50,9	370